

## Econobox™

***A value one and two-bank illuminator designed for alternative care markets.***

PerfectView® is a leading manufacturer of light boxes for medical imaging professionals. The Econobox offers a cost-effective option for the niche industry care provider. Featuring a slim metal wrapper, the Econobox is designed to allow side-by-side mounting with minimal ledge breaks. This viewbox is ideal for private practice locations, offering an economical solution for the cost conscious professional.

### ***Features and Benefits***

- Two 15 watt daylight lamps per viewing area
- UL and CUL Approved
- Instant start, low leakage magnetic ballast
- Two-way rocker on/off switch
- Ribbed ventilation area
- Slim metal wrapper eliminates ledge breaks for side by-side mounting
- Protruded film ledge allows additional viewing space
- Removable plexiglas® panel – no tools required!
- Customized cord mounting (bottom left or right corner)

Also: Ask about our complete line of CPAC Imaging chemistry and silver recovery equipment! PerfectView illuminators are designed by CPAC Imaging in Leicester, NY.




# Specifications

## 1 Bank

Amp: 1.2  
Power rating: 110 VAC power supply  
Dimensions: 20 5/8" H x 14 3/8" W x 3 1/2" D  
Viewing Area: 14 x 17"

## 2 Bank

Amp: 2.4  
Power rating: 110 VAC power supply  
Dimensions: 20 5/8" H x 28 3/4" W x 3 1/2" D  
Viewing Area: (2) 14 x 17"



**PerfectView**<sup>®</sup>  
by CPAC Imaging